

# RHEOSEPT® FD plus

Highly concentrated disinfectant for surfaces and non-invasive medical devices



- economical high concentrate for manual disinfection
- excellent performance at low concentrations
- effective against enveloped, noro-, rota- and adenoviruses
- for higher efficacy requirements
- aldehyde-free

CE 0044 N-113309

## Application areas

Suitable for waterproof surfaces as well as medical furniture, items and equipment. Can be used in bedpan washer-disinfectors. Not suitable for disinfection of skin surfaces.

## Instructions for use

Prepare the working solution from water and surface disinfectant concentrate in accordance with desired concentration and application duration. Apply a thin layer of diluted product to the pre-cleaned surface, ensuring its complete wetting. To remove disinfectant residues from surfaces, wipe them with water after the exposure time. Use freshly diluted solutions only. Only to be used by instructed and trained personnel. Consult the instructions for use of the medical device prior to its disinfecting.

## Efficacy spectrum

Bactericidal, yeasticidal, virucidal (enveloped viruses, adeno-, noro- and rotaviruses).

All properties of the efficacy spectrum have been verified by accredited laboratories.

Efficacy	Norm	Concentration	Contact time
bactericidal/ yeasticidal	EN 13727, EN 13624, EN 16615	0,5/1/3 %	60/15/1 min.
bactericidal/ yeasticidal	VAH	0,5/1/3 %	60/15/5 min.
enveloped viruses	EN 14476	3 %	15 min.
enveloped viruses	EN 14476	1 %	60 min.
enveloped, noro-, adeno-, rotaviruses	EN 14476	5 %	15 min.
norovirus	EN 14476	2 %	5 min.
norovirus	RKI/DVV	1 %	60 min.

## Listings

Listed by the VAH (Verband für Angewandte Hygiene e.V.) and IHO (Industrieverband Hygiene und Oberflächenschutz).

## Registrations

Dual-use item: medical device Class IIa (CE 0044), biocidal product (in Germany, BAuA Reg. Nr. N-113309).

## Composition

100 g solution contains: 5,0 g didecyldimethylammonium-chloride, 1,9 g Benzyl-C12-16-alkyldimethyl, chloride, 4,0 g N-(3-Aminopropyl)-N-dodecylpropan-1,3- diamin. Other components: complexing agents, alcohols, fragrances.

## Physical and chemical properties

pH value: 9 (ready-to-use solution)  
pH value: 11.0-11.5 (concentrate)  
Smell: perfumed  
Colour: yellow  
Form: liquid

## Warning notices /special advice

Causes severe skin burns and eye damage. Very toxic to aquatic organisms. Toxic to aquatic organisms, with long lasting effect.

Please report any serious incidents related to this product to the manufacturer and the competent authority.

Use biocidal products carefully. Always read the label and product information before use.

Store away from light at 5 - 30 °C and closed.

Note the batch number and expiration date printed on the packaging.

## Environmental information

NW-Chemie GmbH produces according to advanced, safe and environmentally friendly processes in compliance with high quality standards.

## Additional information

Further information related to this product can be found at [www.rheosol.com](http://www.rheosol.com). If you have any particular questions, please contact our customer service. Telephone: +49 2241-3923-70, e-mail: [info@rheosol.com](mailto:info@rheosol.com).

## Order information

Reference number (REF)	Item	Content	Sales unit	Pallet quantity
20502352	RHEOSEPT-FD plus bottle	1000 ml	10 pcs/carton	560 pcs
20502050	RHEOSEPT-FD plus canister	5 L	1 pc	56 pcs
20502070	RHEOSEPT-FD plus canister*	5 L	1 pc	69 pcs
20502100	RHEOSEPT-FD plus canister	10 L	1 pc	60 pcs

\*suitable for bedpan washer-disinfectors