

Product information: 20476

bedpan clear rinsing agent

- softening and rinsing agent with good cleaning capabilities in one
- High material compability as well for components with sensitivity for acids (e.g. out of non-ferrous metal).
- no residua on rinsing items and inside the machine (injectors and lines)
- adequate for all thermal bedpan rinsers



RHEOSOL-STN Bedpan clear rinsing agent

Neutral special clear rinsing agent with softener

Product description

Neutral as well as low-foaming rinsing agent and softener for bedpan conditioning in thermal disinfection procedure. For bedpan and faeces rinsers in the medical area. Adequate for thermal conditioning of disinfecting care utensils (high quality steel, plastics, glass) in all commercial machines (appropriate to DIN 15833).

Fields of application

For bedpan and faeces rinsers in the medical field.

Dosing instructions

Please note the directives of the machine manufacturer.

Contact times

Dosing specifications	Dosing details	Remark
recommended	0,5 to 4ml per litre rinsing liquor	depending on machine settings and water hardness

Composition

15-30 % organische Säuren. Weitere Inhaltstoffe: Methylchloroisthiazolinone, Methylisothiazolinone.

Physico-chemical data

- pH-value (application) 7,0 - 7,8

Specific advices

The product is not a dangerous good according to dangerous goods classification. When handling liquid

chemicals the valid preventive measures should strictly be followed. Pay attention to the labelling and the product information prior to use. At least 12 months without particular limitation storable. Store closed, frost-protected and cool.

Forms of delivery

Product code	Package size	Packaging unit
20476050	RHEOSOL-STN Bedpan clear rinsing agent canister 5 litres	canister

Accessories

Product code	Package size	Packaging unit
48017002	Canister tap for 5 and 10-litres-canister with plug 51	piece
48020002	Canister pump plug 51mm, 30 ml pump stroke canister pump for 5 and 10-litres-canister blue or white with UN-code Riser pipe length 300 mm	piece
48038005	Canister key plug 51	piece

Certificate and information

For additional information on this product, please see → www.rheosol.com – For specific questions or information about expert reports, please contact our office staff: Telephone: +49 2241 - 3923-70 e-mail: info@rheosol.de.

Environmental information

Wachendorff-Chemie GmbH produce under application of advanced, safe and environmentally friendly procedures and in compliance with high quality standards.